


Capabilities for high repetition rate streak cameras

Sydor
Instruments



P.A.Jaanimagi, R.Boni and M.Pavia

Sydor Instruments, Rochester NY

SYDOR
INSTRUMENTS, LLC

Single-shot, high precision streak cameras can be operated at 10Hz

*Sydor
Instruments*

- The SNR is reduced at the low end, no impact on the high end
 - fast decay time phosphors are less efficient
- The number of spatial resolution elements is reduced
 - CCD readout is faster with 32x binning in parallel direction
- 20Hz operation possible using a global background frame
 - standard procedure acquires background after each data frame
- Extension to 250Hz with a single channel
 - burst or cw operation

SYDOR
INSTRUMENTS, LLC



Outline


**Sydor
Instruments**

- Review single-shot operation of the ROSS camera.
- What is different for 10Hz operation?
- Options to extend operation to higher frequency.

SYDOR
INSTRUMENTS, LLC



Streak cameras are electron-optical imaging devices that convert the time dimension to a spatial dimension



- A line image is projected onto the photocathode, and the photoelectron replica is imaged onto the screen.
- The polarity of an initial dc bias voltage on the deflection plates is reversed as the electron pulse traverses the region, causing the line image to streak across the screen.

The ROSS camera is a comprehensive diagnostic system with autofocusing and self-calibration capability



High-precision measurements require precise setup and calibration of the diagnostic.

The optical calibration module incorporates an extensive imaging, flat-field correction, geometric-distortion and time-calibration capability.

All functions can be accessed and monitored remotely.

**Rochester Optical Streak System technology
has been transferred to Sydor Instruments**

All pixels in the recorded image were not created equal




The system response is dependent on the spatial and temporal position and the focusing conditions.

To make pretty pictures and FWHM measurements requires minimal calibration.

Serious quantitative measurements, ($\text{SNR} > 10$, $\text{DR} > 10$, nonlinearity $< 10\%$), require extensive calibrations.

These issues are addressed in the ROSS design with the Optical Calibration Module.


The input imaging system is an Offner triplet with motorized controls for the dual object planes and the secondary mirror




Geometric distortion correction of the P820-26 streak tube has spots aligned to 0.49 ± 0.32 pixels



Reticle has $10\mu\text{m}$ openings on $400\mu\text{m}$ centers, LED is pulsed at 5Hz.




Flat-fielding a CCD recording system elevates its performance to a near-quantum-limited detector



Time calibration is accomplished with an on-board 664nm, 2GHz comb generator


Verify system operation prior to shot.

An external source may also be coupled into any or all channels.



4 streak speeds,
(2-25ns), can be
selected remotely


The P820 streak tube provides a 2ps time resolution capability for OMEGA-EP diagnostics




Traps in the CCD can produce a latent image and result in an artificial offset for the background level

Sydor
Instruments

- A fraction (10^{-3} to 10^{-4}) of the CCD electrons are trapped in the wells.
- At -40°C , they are thermally released with a 45 minute time constant.




⇐ Parallel shift direction



The major components of a streak camera must be re-evaluated for high repetition rate operation

*Sydor
Instruments*



Mechanical shutters have limited operational bandwidth

*Sydor
Instruments*

Diameter (mm)	2	3	6	25
Maximum freq. (Hz)	100	50	20	5

Uniblitz LS Series

Irising shutters must be replaced often – disposable item.

Solution:

Use a slot in a rotating disc.

Need sensors and controls to synchronize the opening with the external event.


Photocathodes may become a maintenance item

Sydor
Instruments

Photocathodes are thin, resistive films with limited thermal conductivity.

The main failure mode is loss of QE due to ion feedback from the residual gas in the tube. High rep. rate usage will accelerate the degradation.

Expected lifetime < 5 years.



Standard deflection plate drive circuits also have limited operational bandwidth

Sydor
Instruments

- Avalanche transistor stacks switching $\pm 1\text{kV}$ into 50Ω
 - 10Hz rep. rate OK, at 100Hz the voltage will droop
 - lifetime is unknown, $> 10^7$ cycles

Solution A:


- Semiconductor switches have operated at 1kHz.
 - need $100\mu\text{J}$ of optical energy per switch
 - very low jitter

Solution B:

- Travelling-wave deflection
 - lower voltage and fewer components

Travelling-wave deflection systems have higher bandwidth and higher deflection sensitivity

A 50 Ω stripline in a meander configuration may require only $\pm 100\text{V}$ transients.



Electron beam velocity is in phase with the propagation speed along the stripline.

A fast decay time phosphor and a faster readout CCD must replace the standard ROSS recording system

*Sydor
Instruments*



**P22 has ms decay time +
long low level tail**

P48 has 100ns decay time

Parallel shift overhead:

E2V 42-40BI \Rightarrow 100ms

KAF 1001E \Rightarrow 15ms



The high rep. rate recording system has 5X lower gain

**Sydor
Instruments**

	P22 + E2V 42-40BI	P48 + KAF 1001E
Gain (CCD e⁻ / pe⁻)	100	20
# Pixels	2048 x 2048	1024 x 1024
Pixel size (μm)	13.5	24.0
Array size (mm)	27.65	24.58
Parallel shift (μs)	50	15
N_{read} (e⁻) 400kHz	5	8
I_{dark} (e⁻/pixel/s)	0.05 (-40°C)	0.6 (-30°C)

The streak camera performance at 10Hz is still limited by the streak tube

Sydor
Instruments

Divide the CCD into 2 channels and consider a 20ns streak @
0.1ma / channel \Rightarrow 12500 photoelectrons / channel / 20ps pixel

Binning is 256 x 1 or 2 superpixels / channel.

At 400kHz the 8 electron single pixel read noise increases to
 $N_{\text{read}} \approx 55$ electrons / superpixel with background subtraction.

The rms system noise is equivalent to

4 photoelectrons / channel / 20ps pixel

$\text{SNR}_{\text{peak}} = \text{sqrt}(12500) \approx 100$, dynamic range = $12500 / 4 \approx 3000$

Low end signals are limited by the recording system.

Are there alternative recording systems?

Sydor
Instruments

Frame transfer CCD

- twice the parallel transfers, twice the area

Electron-multiplying CCD

- only small area CCD's available, noise factor

CMOS

- gain non-uniformity, lower QE and higher read noise

Hybrid CCD-CMOS

- new, attractive option

Electron-bombarded CCD

- no phosphor, high gain
- issues with tube processing and HV static discharges
- E.D.Savoye et al., SPIE 203, 59 (1979)

Can we increase the streak camera rep. rate?

*Sydor
Instruments*

Cummulative mode:

- Synchroscan @ 76MHz mode-locked laser frequency
resonant LC circuit to drive the deflection plates

M.C.Adams et al., Opt. Comm. 26, 273 (1978)

- Semiconductor switching at 1kHz, recombination limited
low jitter allowed the detection of ps pulses

G.Mourou and W.Knox, APL 35, 492 (1979)

Circular scan:

- 2 orthogonal sets of deflection plates, RF 90° out of phase
E.K.Zavoiskii and S.D.Fanchenko. Appl. Opt. 4, 1155 (1965)

C.B.Johnson, et al., Appl. Opt. 19, 3491 (1980) – annular Reticon


A CCD camera can be used as a slow scan streak camera during the readout phase. (dual-axis streak)



Record a single channel streak in the first N rows of the CCD

Bin N rows together in the serial register

During digitization, streak the next event into the first N rows



KAF 1001E with N=100

Parallel shift $15\mu\text{s} \times 100 = 1.5\text{ms}$


1000 ADC @ 400kHz = 2.5ms

Single channel streaks @ 250Hz

Operate c.w. or in burst mode collecting a full frame of data

5Hz burst mode operation of a ROSS camera during the CCD readout using the parallel shift to isolate individual streaks

Parallel shift direction ↑



The CCD is being used as a slow sweep streak camera in the parallel direction

Single-shot, high precision streak cameras can be operated at 10Hz

*Sydor
Instruments*

- The SNR is reduced at the low end, no impact on the high end
 - fast decay time phosphors are less efficient
- The number of spatial resolution elements is reduced
 - CCD readout is faster with 32x binning in parallel direction
- 20Hz operation possible using a global background frame
 - standard procedure acquires background after each data frame
- Extension to 250Hz with a single channel
 - burst or cw operation

SYDOR
INSTRUMENTS, LLC

