

UCRL- 92362  
PREPRINT

Cetylpyridinium Chloride - A Versatile  
Reagent for Precipitation Titrations

Walter S. Selig

This paper was prepared for submittal  
to  
J. Chemical Education

March 1985

Lawrence  
Livermore  
National  
Laboratory

This is a preprint of a paper intended for publication in a journal or proceedings. Since changes may be made before publication, this preprint is made available with the understanding that it will not be cited or reproduced without the permission of the author.

VAULT REFERENCE COPY

#### DISCLAIMER

This document was prepared as an account of work sponsored by an agency of the United States Government. Neither the United States Government nor the University of California nor any of their employees, makes any warranty, express or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial products, process, or service by trade name, trademark, manufacturer, or otherwise, does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States Government or the University of California. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States Government or the University of California, and shall not be used for advertising or product endorsement purposes.

## CETYLPYRIDINIUM CHLORIDE - A VERSATILE REAGENT FOR PRECIPITATION TITRATIONS

For the study of potentiometric titrations only very simple equipment is required: a pH/millivolt meter, a reference electrode, and a home-made plastic-coated graphite indicating electrode (1). Quaternary ammonium halides were found very suitable as titrants for precipitation titrations (2). We recommend cetylpyridinium chloride (CPC), a compound which is relatively non-toxic, inexpensive, and very versatile. It can be used for the determination of many inorganic and organic anions. Many cations can be determined after conversion to their halide or cyanide complexes.

### BACKGROUND

Willard and Smith (3) recommended tetraphenylarsonium chloride ( $\phi_4\text{AsCl}$ ) for the gravimetric determination of some large anions and some complex halides. In 1968 Baczuk and DuBois (4) used  $\phi_4\text{AsCl}$  for the potentiometric titration of perchlorate, using a perchlorate ion-selective electrode (ISE). We have found that quaternary ammonium halides can replace  $\phi_4\text{AsCl}$  as the titrant in this and many other titrations (5). In earlier work we used cetyltrimethylammonium bromide as titrant (6). Recent work has shown that CPC is preferable because of its higher solubility in water, which makes it possible to use more concentrated titrant solutions (2). We have also replaced the fairly expensive commercial ISE's with inexpensive home-made plastic-coated graphite indicating electrodes (1).

---

Work performed under the auspices of the U.S. Department of Energy by the Lawrence Livermore National Laboratory under contract No. W-7405-ENG-48.

## EXPERIMENTAL

The preparation of the polyvinyl(chloride)/dioctylphthalate-coated graphite rod was described previously in this Journal (1). The titrant was an aqueous solution of cetylpyridinium chloride which is available from various suppliers.<sup>1</sup> It was prepared by dissolving 3.58 g of the monohydrate in approximately 250 ml of warm water, and diluting to volume with cold distilled water.

The experimental conditions such as feasible pH ranges, sample titration curves, and statistics of recoveries have been described in a review of the analytical applications of quaternary ammonium halides (2) and some recent work (7,8).

A few conditioning runs should precede each experiment. This is also good practice when working with ISE's. In highly acidic solutions the coated graphite indicating electrode will deteriorate more quickly than in other media. This requires re-coating of the graphite sensor as outlined in reference (1).


## INORGANIC ANIONS

A list of inorganic anions that can be determined with CPC is given in Table 1. Many anions can be titrated, with the exception of those of the elements of groups IA through VA and group VIII. Cations in group VIIIA can easily be converted to their halides by adding an excess of potassium bromide or chloride in acid solutions. A large number of complex cyanides can

---

<sup>1</sup>Aldrich catalog no. 85,556-1, Alfa catalog no. 13,651; Eastman catalog no. P5361. The reagent from Alfa was the least expensive material. One liter of a 0.01 M solution costs \$0.17.

similarly be determined, some in acid solutions because of their great stability. A typical reaction, for the determination of  $\text{Pt}^{4+}$ , is given by the equation. In many cases the bromide complex will yield sharper titration



curves than the corresponding chloride because of the lower solubility of the resulting cetylpyridinium precipitate.

Table 1 is not complete and many other applications are no doubt possible. Determinations of the stoichiometry of some of these reactions yielded rather unexpected results. For instance, the titration of  $\text{Bi}^{3+}$  in acid solution in the presence of excess bromide required 1.5 mole of CPC per  $\text{Bi}^{3+}$ . The composition of the resulting precipitate conformed, according to elemental analysis, to the formula  $(\text{C}_{21}\text{H}_{38}\text{N})_3\text{Bi}_2\text{Br}_7$  (8). Molybdates and tungstates form isopoly compounds which are quite sensitive to pH and yield complex compounds in the titration.(9)

#### ORGANIC ANIONS

A list of organic anions that were titratable vs CPC is presented in Table 2. Only water-soluble compounds can be determined because nonaqueous, or partially nonaqueous media will dissolve the organic coating of the graphite indicating electrode. However, many alkali metal salts of organic acids and dyes are water soluble and thus can be determined. It should be noted that a minimum of 10 carbon atoms is required in the analyte for the successful titration of soaps and anionic detergents.

## CONCLUDING REMARKS

We have shown numerous applications of quaternary ammonium halides for the determination of inorganic and organic anions. It is left up to the ingenuity of the student to investigate additional anions and to devise possible separations. Thus:

1.  $\text{Tl}^{3+}$  can be determined as the halide in the presence of  $\text{Tl}^+$ .
2.  $\text{Te}^{4+}$  can be determined as the halide in the presence of  $\text{Te}^{6+}$ .
3.  $\text{Te}^{4+}$  can be determined in the presence of Se.

## LITERATURE CITED

- (1) Selig, W. S., J. Chem. Educ., 61, 80 (1984).
- (2) Selig, W. S., Z. anal. Chem., 312, 419 (1982).
- (3) Willard, H. H., and Smith, G. M., Ind. Eng. Chem. Anal. Ed., 11, 186, 269 (1939).
- (4) Baczuk, R. J., and DuBois, R. J., Anal. Chem., 40, 685 (1968).
- (5) Selig, W. S., Talanta, 26, 1061 (1979).
- (6) Selig, W. S., Mikrochim. Acta, II, 373, 437 (1979).
- (7) Selig, W. S., Mikrochim. Acta, 1984 II, 455.
- (8) Selig, W. S., Z. anal. Chem. (in press).
- (9) Selig, W. S., unpublished data (1982).





Table 1. Some Inorganic Anions Titratable vs CPC, Classified According to Subgroups of the Periodic Table.

VIA	VIIA	VIIIA	IB	IIB	IIIB	IVB	VB	VIB	VIIB	
$\text{Cr}_2\text{O}_7^{2-}$	$\text{MnO}_4^-$	<u>cyanides</u>	<u>halides</u>	$\text{AuCl}_4^-$	$\text{Zn}(\text{CN})_4^{2-}$	$\text{BF}_4^-$	$\text{SnCl}_2^{2-}$	$\text{PF}_6^-$	$\text{S}_2\text{O}_8^{2-}$	$\text{ClO}_4^-$
$\text{MoO}_4^{2-}$	$\text{ReO}_4^-$	$\text{Fe}(\text{CN})_6^{3-}$	$\text{Ru}_2\text{Cl}_{10}^{4-}$	$\text{Au}(\text{CN})_2^-$	$\text{Hg}(\text{CN})_4^{2-}$	$\text{TlCl}_4^-$	$\text{SnCl}_4^{2-}$	$\text{AsF}_6^-$	$\text{TeCl}_6^{2-}$	$\text{IO}_4^-$
$\text{WO}_4^{2-}$		$\text{Fe}(\text{CN})_6^{4-}$	$\text{RhCl}_6^{3-}$		$\text{HgCl}_4^{2-}$	$\text{TlBr}_4^-$		$\text{SbF}_6^-$		$\text{I}^-$
		$\text{Co}(\text{CN})_6^{3-}$	$\text{PdCl}_4^{2-}$		$\text{HgBr}_4^{2-}$			$\text{Bi}_2\text{Br}_7^{3-}$		
		$\text{Ni}(\text{CN})_4^{2-}$	$\text{PdCl}_6^{2-}$							
		$\text{Pd}(\text{CN})_4^{2-}$	$\text{OsCl}_6^{2-}$							
		$\text{Os}(\text{CN})_6^{4-}$	$\text{IrCl}_6^{2-}$							
		$\text{Pt}(\text{CN})_4^{2-}$	$\text{IrCl}_6^{3-}$							
		$\text{Pt}(\text{CN})_6^{2-}$	$\text{PtCl}_4^{2-}$							
			$\text{PtCl}_6^{2-}$							

Table 2 Organic Anions Titratable vs CPC

---

Nitroform

2,4-Dinitrobenzenesulfonate

2,4,5-Trichlorobenzenesulfonate

Picrate

Picrylsulfonate

Tetraphenylborate

Cyanotriphenylborate

Bromophenol blue

Bromocresol purple

Bathophenanthroline disulfonic acid (disodium salt)

Anionic Surfactants and Soaps (sodium salts)

Dioctylsulfosuccinate

Dodecylsulfate

Dodecylbenzenesulfonate

Dodecanesulfonic acid

Octylsulfate

Dodecanoate

Tetradecanoate

Stearate

Oleate

---