

UCRL-92363
PREPRINT

The Sequential Precipitation
Titration of Arsenite and Arsenate
with Silver Nitrate

Walter S. Selig

This paper was prepared for submittal
to
Microchemical Journal

March 1985

Lawrence
Livermore
National
Laboratory

This is a preprint of a paper intended for publication in a journal or proceedings. Since changes may be made before publication, this preprint is made available with the understanding that it will not be cited or reproduced without the permission of the author.

CIRCULATION COPY
SUBJECT TO RECALL
IN TWO WEEKS

DISCLAIMER

This document was prepared as an account of work sponsored by an agency of the United States Government. Neither the United States Government nor the University of California nor any of their employees, makes any warranty, express or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial products, process, or service by trade name, trademark, manufacturer, or otherwise, does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States Government or the University of California. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States Government or the University of California, and shall not be used for advertising or product endorsement purposes.

LETTER TO THE EDITOR

THE SEQUENTIAL PRECIPITATION TITRATION OF ARSENITE AND ARSENATE WITH SILVER NITRATE

Chao and Cheng (1) determined arsenite and arsenate by sequential titration with silver nitrate, using a silver/sulfide ion-selective indicator electrode. The method is practical because the solubility products of the precipitates differ by several orders of magnitude ($-\log K_{sp}$ at 20°C for Ag_3AsO_3 , 31.27 and for Ag_3AsO_4 , 22.19). These authors checked and continually adjusted the pH to 11.0 with dilute sodium hydroxide or nitric acid. In the same paper the authors state that, for the titration of arsenite, a buffer is needed and recommend a sodium hydroxide solution containing some carbonate.

Masson (2) determined arsenate by a similar method using a borate buffer. We have found this buffer (a 0.025 M solution of sodium tetraborate decahydrate) suitable for the sequential microdetermination of arsenite and arsenate and consider the method simpler than the method involving pH adjustment used by Chao and Cheng (1). Arsenite yields a steep titration curve with silver nitrate in aqueous solution. Arsenate, however, because of the larger solubility of its silver salt, yields much shallower titration curves which are enhanced by a partially nonaqueous medium such as 50% methanol. When both anions coexist in solution, the arsenite break is greatly suppressed in 50% methanol. In 25% methanol somewhat better titration curves are obtained for arsenite, but its recovery is not as good as in aqueous

solution, as shown in Table 1. We have at present no explanation for the suppression of the arsenite break in the presence of arsenate in partially nonaqueous media.

A typical titration curve for 0.9 mg of arsenic as arsenite and of 0.6 mg of arsenic as arsenate is shown in Figure 1. Recalculation of the data of Chao and Cheng (1) showed that for arsenite their recoveries varied from -1.02 to +6.1% of the nominal, while for arsenate they varied from -0.9 to +7.6%. This is consistent with our results.

In summary, microamounts of arsenite and arsenate can be determined in borate-buffered aqueous solution by potentiometric titrimetry with silver nitrate.

WALTER S. SELIG

Lawrence Livermore National Laboratory

University of California

Livermore, California 94550

REFERENCES

1. Chao, E. E., and Cheng, K. L., Stepwise titration of some anion mixtures and determination of K_{sp} of silver precipitates with silver ion selective electrode. Anal. Chem. **48**, 267-271 (1976).
2. Masson, M.R., Microdetermination of selenium, tellurium, and arsenic in organic compounds. Mikrochim. Acta **1**, 399-411 (1976).

ACKNOWLEDGEMENTS

This work was performed under the auspices of the U.S. Department of Energy by Lawrence Livermore National Laboratory under contract No. W-7405-Eng-48.




Figure 1. Titration curve for 0.9 mg of As as arsenite and 0.6 mg of As as arsenate (in sequence) in aqueous solution.

TABLE 1

Recovery of Arsenite and Arsenate by Titration with 0.01 N Silver Nitrate

<u>Taken (mmol)</u>		Medium	<u>Mean recovery, percent (standard deviation)</u>		
Arsenite	Arsenate		Arsenite	Arsenate	No. replicates
0.0122	0.008	Aqueous	100.5 (0.67)	98.3 (0.98)	4
0.0122	0.008	25% methanol	103.7 (0.51)	93.5 (1.10)	4