

UCRL-94673
PREPRINT

**Microdetermination of Organic Isocyanates
and Isothiocyanates**

Walter S. Selig

**CIRCULATION COPY
SUBJECT TO RECALL
IN TWO WEEKS**

This paper was prepared for submittal to
Microchemical Journal

June 1986

Lawrence
Livermore
National
Laboratory

This is a preprint of a paper intended for publication in a journal or proceedings. Since changes may be made before publication, this preprint is made available with the understanding that it will not be cited or reproduced without the permission of the author.

DISCLAIMER

This document was prepared as an account of work sponsored by an agency of the United States Government. Neither the United States Government nor the University of California nor any of their employees, makes any warranty, express or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial products, process, or service by trade name, trademark, manufacturer, or otherwise, does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States Government or the University of California. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States Government or the University of California, and shall not be used for advertising or product endorsement purposes.

Microdetermination of Organic Isocyanates and Isothiocyanates*

Walter S. Selig

**Lawrence Livermore National Laboratory
University of California, Livermore, CA 94550**

The method of Vinson for the determination of organic isocyanates and isothiocyanates has been adapted to the microscale. The sample size is approximately 0.05 mequiv. This is particularly useful for high molecular-weight materials.

In the standard method for organic isocyanates and isothiocyanates, the compounds are reacted with n-butylamine in dioxane and the excess amine is then titrated with standard acid (1). This method was adapted to the microscale by Karten and Ma (2). These methods require a reaction time of about 45 min at room temperature. Vinson(3) used a dipolar aprotic solvent, N,N-dimethylformamide (DMF), and thereby shortened the reaction time to 5 min for aromatic compounds and 10 min for aliphatic compounds.

We have applied the latter method on the microscale. Using a 10-ml buret, the following parameters are suggested:

Normality of titrant: 0.01 N HCl

n-butylamine reagent: 150 mg per 100 ml of DMF

Sample size: approximately 0.05 mequiv of isocyanate
(isothiocyanates)

Amount of n-butylamine reagent per sample: 5.00 ml

*Work performed under the auspices of the U.S. Department of Energy by the Lawrence Livermore National Laboratory under contract number W-7405-ENG-48.

Otherwise the procedure is the same as in the macro method (3). We favor, however, a potentiometric endpoint. It should be pointed out that the standardization of the hydrochloric acid in the original procedure is in error: The titrant should be standardized vs tris(hydroxymethyl)-aminomethane and not vs potassium acid phthalate.

Some highly fluorinated aliphatic isocyanates were not soluble in DMF. We dissolved these compounds in 1,1,2-trichloro-1,2,2-trifluorethane. The solutions were compatible with the butylamine/DMF reagent. While on addition of water 2 phases were formed, this did not interfere with the colorimetric endpoint or the electrometric endpoint when the electrodes were immersed in the aqueous phase.

Results for a proprietary fluorinated isocyanate were $4.26 \pm 0.05\%$ NCO for 4 replicates, with a standard deviation of 0.029. The corresponding isocyanate equivalent weight was 985 ± 10 for 4 replicates with a standard deviation of 4. The isocyanate equivalent weight is calculated from % NCO by $(4201.7)/(\% \text{ NCO})$.

REFERENCES

1. Siggia, S., and Hanna, J. G., Determination of organic isocyanates and isothiocyanates. Anal. Chem. 20, 1084 (1984).
2. Karten, B. S., and Ma, T. S., Determination of the isocyanate and isothiocyanate group. Micro and semimicro methods. Microchem. J. 3, 507-514 (1969).
3. Vinson, J. A., Determination of organic isocyanates and isothiocyanates in dimethylformamide. Anal. Chem. 41, 1661-1662 (1969).

0943R