

## ATTOMOLE SENSITIVITY IN ISOTOPE TRACING USING ACCELERATOR MASS SPECTROMETRY.

J.S. Vogel<sup>1</sup> and K.W. TurteltaubBiology and Biotechnology Research Program and <sup>1</sup>Center for Accelerator Mass Spectrometry, Lawrence Livermore National Laboratory, Livermore CA 94551

Accelerator Mass Spectrometry (AMS) facilities at LLNL are being used for ultra-low level tracing of compounds and elements in new applications to life science research. AMS was originally developed as a chronometry tool in earth sciences and archaeology, making special use of small carbon samples for radiocarbon dating. Although <sup>14</sup>C and <sup>3</sup>H are already used for tracing organic compounds, the high efficiency of AMS bestows a very large advantage in throughput and in precision over decay counting (on the order of 10<sup>6</sup>). AMS is unrivaled for detecting <sup>10</sup>Be, <sup>26</sup>Al, <sup>36</sup>Cl, <sup>41</sup>Ca, and <sup>129</sup>I. The lifetimes of these isotopes ensure inefficiency in any method based on decay counting, and traditional mass spectrometry is unable to detect their very low concentrations in natural samples (10<sup>-10</sup> - 10<sup>-15</sup>). The use of AMS in detecting atto- to femto-moles of <sup>14</sup>C-labeled biomolecules from natural systems provides wide dynamic range and direct association of molecular fragments to parent compounds. AMS gives scientists a capability for determining the distributions and fates of chemical compounds and certain elements in natural systems at natural concentrations at a time when such data is required for rational risk determination, drug development, and regulatory justification.

This work was performed under the auspices of the U.S. Dept. of Energy at LLNL under contract no. W-7405-Eng-48.